



Yes, you can.®

➤ *Invacare*® **Action³ Comfort**

A new wheelchair for additional comfort!

The *Invacare*® **Action³ Comfort** now offers optimal lateral positioning and additional seating comfort with its new Dartex® smooth and easily looked after coating. It is the new comfort folding wheelchair.



NEW

New comfort

The *Invacare® Action³ Comfort* offers a 30° reclining backrest fitted as standard, with seat and backrest anatomic preformed soft cushions. In addition, the *Action³ Comfort* now features height and depth adjustable armrests, and legrests with knee protection to provide elderly people with optimal comfort.



Dartex® ergonomic seating

The ergonomic seating (seat and backrest cushions, armrests and headrest can be quickly removed) is covered in Dartex®, a supple and very smooth coating, thus offering unmatched positioning and comfort.

Standard features



Legrests

The multi-adjustable legrests fitted with articulated footplates provide perfect positioning.



Headrest


Multi adjustable; comfort and positioning are optimal



Armrests

The combinations of adjustments ensure ultimate seating positioning.

Technical data




40,5 / 43 / 45,5 / 50,5 cm




45 cm



53* / 50* / 48* cm



56 cm




17 - 26 cm



37 - 49 cm



0° - 30°



59,5 / 62 / 64,5 / 69,5 cm




108 cm




113 cm



86 cm



32 cm



26,5 kg



125 kg

* Seat height including ergonomic cushion (6 cm thickness)

Frame colour



Sapphire Blue



Dark Red



Pearl Grey

Upholstery



Grey Dartex

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited
 South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
 Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016
 E-mail: uk@invacare.com · www.invacare.co.uk
 Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
 Unit 5 · Seatown Business Campus
 Seatown Road · SWORDS · County Dublin
 Tel. (+353) 181 107 84 · Fax (+353) 181 107 85
 Email: ireland@invacare.com · www.invacare.ie