



Yes, you can.[®]

Invacare[®] Ben 9

Simple and practical

The *Invacare[®] Ben 9* is the entry level chair in the Invacare manual range. The lightweight construction of this transit chair makes it easy for even an elderly carer to propel providing a comfortable and adaptable chair.



➤ Added value

Invacare® has built in many beneficial features as standard on the Invacare® **Ben 9**, all designed to make your life easier. These include the hinged backrest and detachable armrests for ease of transfer and to reduce the overall size when folded for transportation. In addition, the footrests are swing away and detachable to aid the transfer process.



➤ Smart and practical

The **Ben 9** is available in a choice of two smart frame colours, with a charcoal black canvas seat and backrest.

A wide range of options are available, from contoured cushions to extended backrests and headrests for added support and comfort. Simply choose the options to suit your individual needs.

Options



Folding backrest

Reduces overall size when folded for transportation or storage.



Swing away footrests


Footplates fold up as shown, and footrests swing away, allowing easier access in or out of the chair.



Foldable cross frame

Lightweight frame folds flat for easier portability.

Technical data




38 - 45 cm




40.6 - 43.2 cm



48 cm



43 cm



22 cm



33 - 43 cm



5°



8°


Frame colour



Black



Red



66 cm




96.5 cm



typ. 91 cm



67 cm



appr. 17 kg



max. 127 kg

Please note colours may vary slightly from those displayed above

Invacare® reserve the right to alter product specification without prior notice. All information believed to be correct at time of print.



INVACARE® Limited · South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
Tel. +44 (0) 1656647327 · Fax +44 (0) 1656649016
E-mail: uk@invacare.com



NHS/Ben 9 - UK - 01/04