



**Yes, you can.<sup>®</sup>**

## **Invacare<sup>®</sup> Duo**

### **Versatile and affordable**

The *Invacare<sup>®</sup> Duo* offers an economical solution to your mobility needs. With its wide range of seat widths up to 24" and 20 stone maximum occupant weight, it makes for a comfortable and adaptable chair.



## Adaptable

This practical wheelchair has been developed from the **Ben 8L** and offers many features as standard to meet client requirements. It's steel construction with 22" quick release wheels and folding backrest makes for a compact and manoeuvrable chair. The **Invacare® Duo** features a rear wheel set back option for increased stability; this can be easily achieved by rotating the rear axle block.



## Practical

The chair is available in two smart frame colours and can be easily adapted to accommodate alternative seating options, depending upon client requirements. A simple cross frame design with swing away footrests and removable armrests allows for easier access in or out of the chair. For ease of storage and transportation the chair will dismantle and fold away.

### Features



22" quick release solid wheels




Removable armrests



Easier side transfers

### Technical data




40 - 62.5 cm




42 - 47 cm



48 cm



47 cm



22 cm



35 - 44 cm



3.5°



10°


#### Frame colours



Red



Black



63 - 85 cm




68 - 72 cm (Folded)



104 cm



78 cm



19.3 - 19.7 kg



max. 127 kg

Invacare® reserve the right to alter product specification without prior notice. All information believed to be correct at time of print.



INVACARE® Limited · South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY  
Tel. +44 (0) 1656647327 · Fax +44 (0) 1656649016  
E-mail: uk@invacare.com

