



Yes, you can.®

## ➤ Invacare® **Spectra Blitz**

### **Go and grow with the Invacare® Spectra Blitz**

The *Invacare® Spectra Blitz* is a versatile power chair for children. Its compact design provides excellent mobility and driving performance. The seat is highly adaptable and allows precise adjustment for optimal comfort and posture management.





Height adjustable back rest.



Adjustable seat width from 25 to 39 cm.



Settable back rest angle for easy positioning.



Rear wheel suspension for better shock absorption.



Angle adjustable leg rests which neatly swing away.




Easy control with the new « Shark II » electronics.

### Technical data

25 - 39 cm	25 - 38 cm	45 cm	40 - 47 cm	-5° - 25°
56 cm	94 cm	91 cm	abt. 62 cm	5 cm
16 %	2 x 40 Ah	2 x 106 W	up to 21 km	6 km/h
79,5 kg (incl. batteries)	75 kg			

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

### Frame colours



INVACARE® Limited  
 South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY  
 Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016  
 E-mail: uk@invacare.com · [www.invacare.co.uk](http://www.invacare.co.uk)  
 Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd  
 Unit 5 · Seatown Business Campus  
 Seatown Road · SWORDS · County Dublin  
 Tel. (+353) 181 107 84 · Fax (+353) 181 107 85  
 Email: ireland@invacare.com · [www.invacare.ie](http://www.invacare.ie)

Spectra Blitz - UK - 07/07 - Art. No. 1450089  
 © Invacare® International SARL 2007. All rights reserved.

