

## Technical data

Standard Seat	44 - 48 cm	41 / 46 / 51 cm	45 cm * <sup>1</sup>	48 / 54 cm * <sup>2</sup>	23 - 32 cm	29 - 47 cm * <sup>3</sup>	0 - 50°
Contour Seat	44 - 53 cm	41 / 46 / 51 cm	45 cm * <sup>1</sup>	64 cm * <sup>2</sup>	23 - 32 cm	29 - 47 cm * <sup>3</sup>	0 - 50°
Flex Seat	44 - 53 cm	38 / 43 / 48 cm	45 cm * <sup>1</sup>	55 cm * <sup>2</sup>	23 - 32 cm	29 - 47 cm * <sup>3</sup>	0 - 50°
Standard Seat	0 - 35°	64 cm	98 cm * <sup>4</sup>	98 cm * <sup>5</sup>	75 cm	138 / 165 kg * <sup>6</sup>	180 kg
Contour Seat	0 - 35°	64 cm	108 cm * <sup>4</sup>	98 cm * <sup>5</sup>	75 cm	138 / 165 kg * <sup>6</sup>	180 kg
Flex Seat	0 - 35°	64 cm	108 cm * <sup>4</sup>	98 cm * <sup>5</sup>	75 cm	138 / 165 kg * <sup>6</sup>	180 kg
	73,6 Ah (C20)	56 cm	7,5 cm	9° / 16 %	2 x 340 W	39 km	8 km/h

- \*1 ground to seat plate, w/o seat cushion w/o CG tilt
- \*2 measured from seat plate
- \*3 measured from foot plate to seat plate
- \*4 standard configuration, varies with options
- \*5 leg rests at largest angle
- \*6 standard configuration, varies with options w/o CG tilt

All measurements have been carried out under Invacare® standard test conditions, unless otherwise stated:  
 User weight: 100 kg, battery: 73 Ah, speed version: 8 km/h, std seat, seat tilt angle: 4°, seat height: 45 cm, back rest angle 10°

### Frame colours



Electric Red

All information quoted is believed to be correct at time of print.  
 Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited  
 South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY  
 Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016  
 E-mail: uk@invacare.com · [www.invacare.co.uk](http://www.invacare.co.uk)  
 Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd  
 Unit 5 · Seatown Business Campus  
 Seatown Road · SWORDS · County Dublin  
 Tel. (+353) 181 107 84 · Fax (+353) 181 107 85  
 Email: ireland@invacare.com · [www.invacare.ie](http://www.invacare.ie)





**Yes, you can.<sup>®</sup>**

**Invacare<sup>®</sup> TDX<sup>®</sup> SP**

**Get more than mobility**

The new *Invacare<sup>®</sup> TDX<sup>®</sup> SP* is a high-performance power chair for those that want to get the maximum out of their wheelchair. The *TDX<sup>®</sup> SP* features Centre Wheel Drive technology for outstanding manoeuvrability. Active users will appreciate the Stability Lock function which offers maximum stability and safety on uneven terrain and SureStep<sup>®</sup>, for smooth travel over obstacles. Our pioneering tilt-in-space function provides optimum pressure relief and positioning. Even heavier users can enjoy speedy, secure mobility thanks to a powerful drive mode.





### ➤ **Yes to mobility, no to compromise**

The tightest possible turning radius is based in the Centre Wheel Drive construction. This type of wheel base provides excellent manoeuvrability and an intuitive driving experience. The *Invacare*® **TDX**® **SP** is the ideal companion for both indoor and outdoor conditions thanks to good stability and high traction on uneven ground.

### ➤ **Tried, tested .... and tough**

For maximum stability and safety, the patented Stability Lock function keeps all six wheels on the ground. The SureStep® feature enables a smooth ride over even the toughest of obstacles. For transportation, the crash-tested 'tie-down' ensures a secure transportation.

### ➤ **Sit back ... and relax**

Each one of our many seat choices can be set to meet individual requirements. Both seat and backrest angle can be adjusted to suit driving or relaxing. The tilt-in-space function, with its angle of more than 50°, offers optimum positions for both comfort and pressure relief.

### ➤ **Designed for heavier users**

Thanks to sophisticated engineering, robust construction and powerful motors, the **TDX**® **SP** is ideal for heavy-duty use and can support a weight of up to 180 kg.



## Features & Options



### Centre Wheel Drive

This specific driving technology allows excellent manoeuvrability in narrow spaces and offers a responsive driving experience.



### SureStep®

Original patented technology carries the wheelchair up and over obstacles up to 7.5 cm high, and guarantees a smooth downwards transition.



### Stability Lock

For maximum stability and safety, this innovative feature keeps all six wheels on the ground and also has a shock-absorbing effect.



### Tilt-in-space

A maximum seat angle of 50° ensures a comfortable sitting position and aids in pressure relieving.



### Reclining back rest

A maximum backrest angle of 35° provides excellent seating and back support.



### 4-point tie-down

Four strong hooks ensure a stable position during transportation. The **TDX® SP** has been successfully crash-tested.