



Yes, you can.[®]

Invacare[®] Zipper 2

Comfort and reliability

The strong but lightweight construction and host of options and accessories make the *Invacare[®] Zipper 2* easy to configure to suit most requirements. Available in both the transit and self-propel versions, this competitively priced chair is the ideal choice.





▶ Value for money

The *Invacare*® **Zipper 2** has been priced to be highly competitive with similar models on the market, offering value for money together with built-in quality. The **Zipper 2** has been designed specifically with economy in mind, helping to ease pressure on ever tightening budgets, while still giving the flexibility and quality expected by Wheelchair Centre customers. It is also ideal for the voucher scheme.

▶ Compatibility and choice

The **Zipper 2** has been designed around standard NHS specifications, making it is easy to reconfigure and maintain using the many readily-available options and accessories. The ability to take the NHS standard tray when fitted with the full length arm illustrate the flexibility and compatibility of the **Zipper 2**. There are six seat widths and four seat depths to choose from, making the **Zipper 2** the most configurable chair in its class.

▶ Lightweight and portable

Using lightweight materials, the **Zipper 2** weighs less than many other manual wheelchair and folds flat, making it easy to stow in the car, or take with you wherever you are going. The lightweight construction has not, however, compromised strength. The one-piece frame design gives a sturdy and robust chair, whilst offering a capacity of 114kg.

▶ Safety and security

As with all our chairs, safety when driving or when in a vehicle is of paramount importance to us. The **Zipper 2** has been subjected to strict testing to comply with the latest requirements applicable to the transportation of wheelchair occupants within vehicles. Both the transit and self propelling models have been crash-tested using a 4 point tie down system, and have successfully passed all tests. This means that users can have complete confidence in the **Zipper 2**, and in its abilities as a safe, practical lightweight chair.



Options



Head rest extension

Offering extra support to the head and neck.



Transit amputee bracket

For use with transit rear wheels.



Self propel amputee bracket

For use with standard rear wheels.



Armrests

Range of armrest styles for comfort and support.



Transit wheel options

For use with an assistant or carer.



Rear wheel options

For self propelling mode.



Anti-tippers

Giving added stability when travelling on uneven ground.



Tipping lever

For use when kerb climbing.




Calf strap

Gives support and added comfort to the lower legs.


Technical data



Self Propel	35.5 - 48 cm	39.5 - 47 cm	50 cm	42 cm	23 cm	38 - 49 cm
Transit	35.5 - 48 cm	39.5 - 47 cm	50 cm	42 cm	23 cm	38 - 49 cm



Self Propel	5°	7° bend	53 - 65,5 cm	95 cm	108 cm	82,5 cm
Transit	5°	7° bend	48 - 60,5 cm	95 cm	97,5 cm	72 cm



Given values may vary with configurations

Self Propel	16,8 kg	114 kg
Transit	14,5 kg	114 kg

Frame colour



Green



Blue



Celt Blue



Black

Please notice colors may vary slightly from those displayed above.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.