


# GH 2 TRACK HOIST SYSTEM



- Guldman Hoist Systems
- Site surveys
- Sling Assessments

**Kingkraft**  
KINGKRAFT

# Guldmann®



The Guldmann system is a modular system, easy to handle during transportation and installation – and suitable for recycling.

The modules can be combined to form small or larger track systems, room covering systems and combinations of these.

The rails can be mounted on the ceiling and on the wall, or can be free standing and moveable systems.

The hoist can be moved smoothly and easily from one rail component to the other.

Mounting brackets are available for any possible situation and travelling trolleys run quietly in the rail system using very little power.



# Kingkraft

*Please contact us for further information*

Kingkraft Limited, 26D Orgreave Crescent, Dore House Industrial Estate, Sheffield S13 9NQ  
Telephone: 0114 269 0697 Fax: 0114 269 5145 Email: [info@kingkraft.co.uk](mailto:info@kingkraft.co.uk) Web: [www.kingkraft.co.uk](http://www.kingkraft.co.uk)

Due to the process of continuous development we reserve the right to alter specification without notice.