

Uniclass L72101:N5352	EPIC H122:Y311
CI/SfB (74.22)	X (U35)
MB/10/2007	

# MULTIBATH

*Variable height bathing & showering*




The Multibath from Kingkraft combines conventional bath looks with generous dimensions to provide a more complete assisted-bathing system with multiple applications. The internal dimensions of the tub allow for the optional Kingkraft cushioning and bather support system to be used for the comfort of the bather.

Benefits to the user are:

- Height adjustable to suit Carer
- Generous internal shape and size – depth 400mm
- Suitable for use with ceiling and most mobile hoists
- Two model lengths
- Designed to be used with the Kingkraft range of baths supports and cushions
- Flush fit safety fill spout
- Recessed side panel for improved carer position

The styling of the sides of the Multibath means that a conventional bath height is achieved at the lowest level to allow for more independent access by ambulant users. The bath can thus be used by a varied group of users from low to high dependency with either a mobile or ceiling hoist as a mode of transfer.



Model No.	A Length	B Width	C Height approx	Internal length	Internal width
1700	1660	780	530-930	1400	590 max.
2000	2020	780	530-930	1700	590 max.

All dimensions expressed in millimetres



**Maximum Load** : 130kg/21st

Also available is the **Multibath-Lite 1700 Model**, which uses less water for rapid fill and emptying.  
 Depth 280mm  
 Internal Length 1400mm max



Kingkraft Limited, 26D Orgreave Crescent, Dore House Industrial Estate, Sheffield S13 9NQ  
 Tel: 0114 269 0697 Fax: 0114 269 5145 Email: info@kingkraft.co.uk www.kingkraft.co.uk

Due to the process of continuous development we reserve the right to alter specification without notice.



**standard features:**

- Raise and lower
- Integral shower
- Thermostatic mixer for bath fill and shower

**optional extras:**

- Twin waste
- Whirlpool system
- Kingkraft support and positioning systems
- Non standard colours
- Wall mounted mixer valves for variable temperature control
- Different styles and configurations of taps
- Changing table top (either hinged or removable)