

VARIABLE HEIGHT BASIN



The concept of the Kingkraft Variable Height Basin is simple: the standard basin comes complete with mirror, light and mixer tap, with a wide basin surround to provide plentiful arm support where required. However the basin can be designed by the user and carer to provide an unlimited combination of styles.

- Selection of basin shapes, sizes and colours custom designed to suit the individual needs of the user
- Slim basin surround
- Robust stainless steel frame construction
- Variety of tap options eg. Infra Red

*Variable height – 400mm travel*

