

# ROPOX NURSING TABLE

*Change, wash and bathe all on one table*



FAMILY  
**Kingkraft**



A wall mounted nursing table with sink and furnished with two wire baskets underneath the tabletop. The nursing bench is mounted on a wall console. The height is indefinitely adjustable from 30-100cm by means of a hand control in a spiral cord. The nursing table with sink is delivered with cover plate, mixer tap with scalding protection, flexible connecting hoses and draining hose with water trap. A range of sizes available.

**Accessories include:**

- Mattresses
- Bed guard

**Height:** 30-100 cm

**Max. load:** 100 kg

**Sink:** Stainless steel, 40x34 cm

**Models:**

- 40-25591 size 90x73 cm
- 40-25592 size 90x73 cm
- 40-25593 size 140x73 cm w/sink
- 40-25595 size 90x80 cm
- 40-25596 size 90x80 cm
- 40-25597 size 130x80 cm w/sink
- 40-25598 size 130x80 cm

**Wall console:** Height 114 cm/Width 71 cm  
Distance from wall to tabletop: 10 cm

**Mattresses:**

- 40-25793 size 90x70x3 cm
- 40-25792 size 90x75x3 cm
- 40-25797 size 80x75x3 cm

**Safety device:**

Safety stop mounted underneath the table top

**Side covering:**

- 40-25785 95x87 cm for tables, depth 80 cm
- 40-25786 95x80 cm for tables, depth 73 cm

Recommended to avoid squeezing when more nursery tables are mounted placed beside each other

**Bed guard:**

- 40-25960 Foldable, 110 cm  
For nursing table with sink
- 40-25965 Foldable, 60 cm  
For nursing table without sink

