

# LIFT CHAIRS



Wing Back  
Lift Chair

Interchangeable  
Accent Panels



## FEATURES

- Integrated battery backup in case of mains failure
- Heavy-duty motor/actuator
- Up to 113 kg (18 STONE) weight capacity
- Quiet & smooth lift and recline system
- Wing back design offers head support during sleep
- Removable fabric panel on front of arms can be interchanged with any of 3 wooden accents in Oak, Pine or Mahogany (wooden accents included with lift chair)
- Fabrics tested to Gigorette (0), Match (1) and Crib (5)
- Two year limited warranty

## SPECIFICATIONS

Lift Chair Model	Back Width	Seat to Floor	Seat Depth	Seat Width	Top of Back to Seat
C1 Petite	51 cm (20")	41 cm-46 cm (16"-18")	46 cm-51 cm (18"-20")	48 cm (19")	66 cm (26")
C15 Cameo	71 cm (28")	48 cm (19")	48 cm (19")	53 cm (21")	68.5 cm (27")
D30 Comfort	76 cm (30")	46 cm (18")	53 cm (21")	55 cm (22")	66 cm (26")
Wing Back	48.25 cm (19")	42 cm (16.5")	49.5 cm (19.5")	47 cm (18.5")	67 cm (26.5")
560 Mini Lounger	71 cm (28")	48 cm (19")	48 cm (19")	48 cm (19")	63.5 cm (25")
C5 Steel Chair	71 cm (28")	46 cm-51 cm 18"-20"	46 cm-51 cm 18"-20"	56 cm (22")	66 cm-71 cm 26"-28"
C6 Steel Chair	71 cm (28")	48 cm (19")	53 cm (21")	56 cm (22")	66 cm-71 cm 26"-28"
660 Mini Lounger Duet	71 cm (25")	48 cm (19")	51 cm (20")	48 cm (19")	63.5 cm (25")
670 Chairbed	76 cm (30")	51 cm (20")	53 cm (21")	51 cm (20")	68.5 cm (27")



1) Only available on the C15, D30 and Wing Back

"The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice."

UK-LiftChair-Range-6-05-08

**Pride**  
Mobility Products Ltd.  
www.pridemobility.com

# LIFT CHAIRS

THE ULTIMATE IN COMFORT AND STYLE®



Place dealer stamp here



www.pridemobility.com



# LIFT CHAIRS

# LIFT CHAIR FEATURES

# LIFT CHAIRS

C-1 - Petite



- C-1 - Petite**
- Quiet and smooth 3-position lift system
  - All steel frame construction
  - Weight capacity of 113 kg (18 Stone)
  - Compact dimensions make it perfect for the petite individual
  - Adjustable seat-to-floor height
  - Adjustable seat depth
  - All upholstery is removable

- C-5 - Steel Chair**
- 3-way action
  - All steel frame construction
  - 170 kg (27 Stone) weight capacity
  - Adjustable seat-to-floor height and seat depth
  - All upholstery components are fully removable

- C-6 - Steel Chair**
- Infinite number of recline & lounging positions
  - All steel frame construction
  - 136 kg (21 Stone) weight capacity
  - All upholstery components are fully removable

- C-15 - Cameo**
- Quiet and smooth 3-position lift system
  - Comfortable blown-fibre, button-back design
  - Weight capacity of 147 kg (23 Stone)
  - Stylish design that adds to any decor

- D30 - Comfort**
- 3-way action
  - Engineered Furniture Laminate/Hardwood frame
  - 170 kg (27 Stone) weight capacity
  - Head and arm covers

- 560 - Mini Lounger**
- Removable seat and back cushions
  - Engineered Furniture Laminate/Hardwood frame
  - 170 kg (27 Stone) weight capacity
  - Head and arm covers

- 660 - Mini Lounger Duet**
- Infinite number of recline & lounging positions
  - Removable seat and back cushions
  - Engineered Furniture Laminate/Hardwood frame
  - 170 kg (27 Stone) weight capacity
  - Two independent motors
  - Head and arm covers

- 670 - Chairbed**
- Infinite number of recline & lounging positions
  - Removable seat and back cushions
  - Engineered Furniture Laminate/Hardwood frame
  - 170 kg (27 Stone) weight capacity
  - Two independent motors
  - Head and arm covers

Note: There is a 2-year limited warranty on all Pride® Lift Chairs

D30 - Comfort



560 - Mini Lounger



660 - Mini Lounger Duet



670 - Chairbed



C-5 - Steel Chair

C-6 - Steel Chair



C-15 - Cameo

