

Pride® Lift Chairs
The Ultimate in Comfort & Style

LC-101

Lift Chair



A 3-Position Lift Chair that delivers Superb Value and Excellent Comfort

Pride's LC-101 delivers superb value, excellent comfort and exceptional quality. The LC-101 features quiet and smooth operation and stylish designs to add to any home décor.

SPECIFICATIONS

LC-101

Seat to Floor:	47 cm (18.5")
Seat Depth:	51 cm (20")
Seat Width:	51 cm (20")
Top of Back to Seat:	70 cm (27.5")
Distance From Wall:	45.75 cm (18")
Overall Chair Width:	83 cm (32.75")
Weight Capacity:	147 kg (23 Stone)
Back Style:	T-back

FABRIC SELECTION



Vino Pine Blue Almond Oatmeal

Note: Measurements may vary +/- 2.5 cm (1").

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

The Ultimate in Comfort & Style

BUILT WITH PRIDE

Pride
Mobility Products Ltd.

32 Wedgwood Road
Bicester, Oxfordshire OX26 4UL

Tel: 01869 324 600
Fax: 01869 323 070
www.pride-mobility.co.uk

