

# Jazzy® GO-2

The Jazzy® Go-2 delivers high performance manoeuvrability and efficiency in a compact, stylish package. A powerful PG GC programmable controller, in-line motors, and enhanced battery size are just a few of the features that make the Jazzy Go-2 an excellent overall value.



Jazzy® Go-2  
shown in Red

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.com](http://www.pridemobility.com)

BUILT WITH PRIDE  
**Pride**  
Mobility Products Ltd.

## SPECIFICATIONS

<b>Weight capacity</b>	136 kg (21 Stone)
<b>Maximum speed</b>	Up to 5.5 km/h (3.5 mph)
<b>Ground clearance</b>	4.45 cm (1.75")
<b>Turning radius</b>	58 cm (22.75")
<b>Overall length</b>	96 cm (37.75")
<b>Overall width</b>	55 cm (21.5")
<b>Wheels</b>	Drive: 20 cm (8") solid
<b>Front anti-tips</b>	13 cm (5") solid
<b>Rear casters</b>	8 cm (3") solid
<b>Suspension</b>	Limited
<b>Drive train</b>	Two-motor, in-line, rear-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC controller
<b>Specialty controls</b>	No
<b>Battery charger</b>	On-board, 4 amp (standard)
<b>Per-charge range</b>	Up to 24 km (Up to 15 miles)
<b>Battery requirements</b>	(2) 12 volt, 30 AH
<b>Battery size</b>	30 AH
<b>Battery weight</b>	11 kg (24.5 lbs.) each
<b>Base weight</b>	32 kg (70 lbs.)
<b>Standard seat weight</b>	14 kg (30 lbs.)
<b>Max. seat dimensions</b>	46 cm x 41 cm (18" x 16")
<b>Warranty</b>	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

## COLOR CHOICES



Red



Blue

## FEATURES

- 30 AH batteries for greater per charge range
- In-line motor technology for increased torque and efficiency
- Tight turning radius of 58 cm (22.75") for tight space manoeuvrability
- Solid tires for worry free ride
- 40 amp PG GC controller
- 136 kg (21 Stone) weight capacity
- 55 cm (21.5") overall width
- 5 amp off-board charger
- Full length padded armrests

## OPTIONS

- Weather cover
- Cup holder
- Rear basket

**Jazzy® Go-2**  
shown in Blue



**Pride**<sup>®</sup>  
Mobility Products Ltd.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; due to the manufacturers commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.