

# GO-CHAIR®

## TRAVEL CHAIR

The Ultimate Travel Chair from the first name in travel mobility, the Go-Chair® is the first travel chair to combine super portability, excellent manoeuvrability and style at a low, low price. The Go-Chair easily disassembles for transport and fits in most car boots.



Go-Chair®  
Travel Chair  
shown in Red

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.com](http://www.pridemobility.com)

BUILT WITH PRIDE  
**Pride**  
Mobility Products Ltd.

## SPECIFICATIONS

<b>Weight capacity</b>	113 kg (17 Stone)
<b>Maximum speed</b>	Up to 5.5 km/h (3.5 mph)
<b>Ground clearance</b>	4 cm (1.75")
<b>Turning radius</b>	57 cm (22.5")
<b>Overall length</b>	76 cm (30")
<b>Overall width</b>	47 cm (18.5")
<b>Wheels</b>	Drive: 20 cm (8") solid
<b>Front castors</b>	13 cm (5") solid
<b>Rear anti-tips</b>	8 cm (3") solid
<b>Suspension</b>	Limited
<b>Drive train</b>	Two-motor, rear-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC2 controller (non-programmable)
<b>Specialty controls</b>	No
<b>Battery charger</b>	Off-board, 2 amp (standard)
<b>Per-charge range</b>	Up to 13 km (Up to 8 miles)
<b>Battery requirements</b>	(2) 12 volt
<b>Battery size</b>	12 AH
<b>Battery weight</b>	9 kg (20 lbs.) battery pack
<b>Base weight</b>	27 kg (59.5 lbs.) front & rear assembled w/o batteries
<b>Standard seat weight</b>	13.5 kg (30 lbs.) foldable, compact, lightweight seat
<b>Max. seat dimensions</b>	46 cm x 46 cm (18" x 18")
<b>Warranty</b>	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

## COLOR CHOICES



Red



Blue

## FEATURES

- Easy to remove battery pack and compact, lightweight frame
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle boots
- Top speed of 5.5 km/h (3.5 mph) and enough power to get you where you want to go
- Twin motors for traction and precision handling
- Tight turning radius of 57 cm (22.5")



**Go-Chair®  
Travel Chair**  
shown in Blue



**Pride®**  
Mobility Products Ltd.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL

(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

