

Jazzy®

1121

Jazzy®

1121



JAZZY® 1121
Shown in Blue with
Euro-Comfort Seat,
swing-away leg rests
and PG Remote Plus
Controller

The Jazzy® 1121 with enhanced Active-Trac® Suspension (ATS) makes easy work of rough and uneven outdoor terrain. Plus, an articulating rear castor beam and patented resilient, adjustable front anti-tips further serve to add to the 1121's comfort and stability. Even tight-space indoor manoeuvrability is made easy with the 1121's compact overall size and mid-wheel drive design. The Jazzy 1121 with ATS is the perfect power chair for your diverse and demanding lifestyle.



JAZZY® 1121
Shown in Red with
Economy Seat, footboard
and VSI controller

Pride
Mobility Products Ltd.
www.pridemobility.com

32 Wedgwood Road
Bicester, Oxfordshire, OX26 4UL
Tel: 01869 324 600
Fax: 01869 323 070



Specifications:

WEIGHT CAPACITY:	136kg (21 STONE)
MAXIMUM SPEED:	Up to 6 km/h (4 mph)
GROUND CLEARANCE:	9cm (3.5")
TURNING RADIUS:	57cm (22.5")
LENGTH:	97cm (38.125")
WIDTH:	55cm (21.5")
DRIVE WHEELS:	35.5cm (14") pneumatic
CASTORS:	20cm (8") solid
ANTI-TIPS:	15cm (6") solid
SUSPENSION:	Active-Trac
DRIVE TRAIN:	2-motor, mid-wheel drive
BRAKING SYSTEM:	Intelligent braking (electronic, regenerative disc brakes)
ELECTRONICS:	50 Amp VSI
STANDARD SEATING:	High Back with headrest, grey vinyl
STANDARD FOOTREST:	Footboard
BATTERY CHARGER	Off-board, 8 amp (standard)
PER-CHARGE RANGE:**	Up to 40km (25 miles)
BATTERIES:	Two 12 volt 55 amp
BATTERY WEIGHT:	18kg (38lbs.) each
BASE WEIGHT:	61kg (135.5 lbs.)
SEAT WEIGHT:	22kg (49lbs.)
MAX. SEATING	
DIMENSIONS:	46cm x 46cm (18"x18") (Economy)
WARRANTY:	Lifetime limited warranty on frame 1-year warranty on electronics 18-month warranty on drive motor



■ VSI integral controller



■ Easy accessible side charger port



■ Rear articulating castor beam

Features:



- Active-Trac Suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- Rear articulating castor beam for great performance and durability
- One-piece, 5 cm (2") diameter tubular frame for maximum strength and enhanced performance
- Easy front access to batteries
- Easy accessible side charger port
- VSI integral controller

Options:

- Euro Comfort seat
- Solid seat pan
- Solid tyres
- Cup holder
- Pilot Plus controller
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Lap belt
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Power elevating seat
- Rear basket
- Angle-adjustable footplates
- 70° and 80° heavy-duty drop-in footrests

Colours:



Candy Apple Red (std)



Viper Blue (std)



Forest Green (option)

■ Swing-away joystick



■ Elevating leg rests



■ Crutch holder



**The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.*

***Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.*