

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 provides maximum front stability, combined with the advantages of in-line motor technology. Active-Trac[®] Suspension and 6" front and rear casters offer superior performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.



Jazzy[®] Select[™] 6
shown in Red

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.com

BUILT WITH PRIDE
Pride[®]
Mobility Products Ltd.

SPECIFICATIONS

Weight capacity	136 kg (21 Stone)
Maximum speed	Up to 6 km/h (4 mph)
Ground clearance	5 cm (2")
Turning radius	52 cm (20.5")
Overall length	90.5 cm (35.625")
Overall width	60 cm (23.5")
Wheels	Drive: 25 cm (10") solid
Anti-tips	N/A
Castors	15 cm (6") solid (front and rear)
Suspension	Active-Trac®
Drive train	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC/ 60 amp VR2*
Specialty controls	No
Battery charger	Off-board, 5 amp (standard)
Per-charge range	Up to 19 km (Up to 12 miles)
Battery requirements	(2) 12 volt deep cycle
Battery size	30 AH
Battery weight	11 kg (24.5 lbs.) each
Base weight	51 kg (112 lbs.)
Standard seat weight	19.25 kg (42.5 lbs.) high back
Max. seat dimensions	46 cm x 46 cm (18" x 18")
Warranty	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

COLOR CHOICES



Red



Blue

* VR2 only available with power seat option

FEATURES

- 15 cm (6") castors on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- 40 amp, PG GC controller/ 60 amp VR2 (some models)
- Articulating rear castor beam
- Battery removal through side panel
- Easy access charger port and circuit breaker
- Service friendly access to controller
- Height and depth-adjustable foot platform

OPTIONS

- Power elevating seat
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

Jazzy® Select™ 6

shown in Red with power elevating seat option



Pride
Mobility Products Ltd.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL

(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.com

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

