

# MAXIMA

**NOW AT  
35.7 STONE  
WEIGHT CAPACITY**



Maxima 3-wheel scooter  
Shown in Red

The Maxima Scooter combines outstanding performance and elegant styling together with rugged construction for an increased weight capacity of up to 227 kg (35.7 Stone)! With a list of standard features like an easily adjustable tiller, fully directional light package, 56 cm (22") wide medium-back seat, 7-piece disassembly, front anti-tip wheels, ultra heavy-duty drive train, 55 Amp battery capacity and enhanced rear suspension—the Maxima provides heavy-duty brawn with all of the convenience and comfort of a luxury scooter.



Maxima 4-wheel scooter  
Shown in Blue

**Pride**  
Mobility Products Ltd.

Tel: 01869 233 100  
Fax: 01869 233 400



[www.pridemobility.com](http://www.pridemobility.com)

## Specifications:

|                                      |  |
|--------------------------------------|--|
| WEIGHT CAPACITY:                     | 227 kg (35.7 STONE)                                |
| MAXIMUM SPEED:                       | Up to 6 km/h (4 mph)                               |
| GROUND CLEARANCE:                    | 11 cm (4.5")                                       |
| TURNING RADIUS:                      | 3 wh: 132 cm (52")<br>4 wh: 137 cm (54")           |
| LENGTH:                              | 3 wh: 117.5 cm (46.25")<br>4 wh: 119 cm (46.875")  |
| WIDTH:                               | 65 cm (25.5")                                      |
| FRONT TYRES:                         | 8 cm x 25 cm (3" x 10" solid)                      |
| REAR TYRES:                          | 10 cm x 25 cm (4" x 10" solid)                     |
| ANTI-TIPS:                           | 3 wh: Front & Rear<br>4 wh: Rear                   |
| STANDARD SEAT:                       | Compact foldable<br>moulded plastic                |
| SUSPENSION:                          | Rear   |
| BATTERY CHARGER:                     | Off-board, 5 amp                                   |
| PER-CHARGE RANGE:**                  | Up to 32 km (20 miles)                             |
| BATTERIES:                           | (2) 12V deep-cycle, 55 AH                          |
| BATTERY WEIGHT:                      | 17 kg (38 lbs.) each                               |
| TOTAL WEIGHT w/o<br>BATTERIES:       | 3 wh: 79 kg (175 lbs.)<br>4 wh: 83.5 kg (184 lbs.) |
| HEAVIEST PIECE<br>WHEN DISASSEMBLED: | 33 kg (72.5 lbs.) rear section                     |

## Colours:



Red



Blue

## Features:

- 227 kg (35.7 Stone) weight capacity
- Full directional light package
- Ultra heavy-duty drive train with 110-amp controller
- 65 cm (25.5") width for optimal stability
- Flat-free foam rubber-filled tires
- Maximum speed up to 6 km/h (4 mph)
- Front basket
- Ergonomically enhanced 56 cm (22") wide x 47 cm (18.5") deep seating
- Durable and stylish rear polyurethane bumper
- User-friendly console
- Stylish micro-mag alloy wheels
- Standard Pride® warranty

## Options:

- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder

\*The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.\*

\*\*Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.\*