



BOSTON HOME CARE BED

The Boston range of beds can be specified in two different combinations for maximum flexibility and safety. Simply choose either a bed with four plastic feet or with castors.



KEY FEATURES:

- Maximum user weight 180kg (28 stone)
- Specifically designed, static & dynamic mattresses available (See Surfaces section - Page 40)
- Full range of bed accessories (See Page 36)

ACCESSORIES:

IV (Transfusion Pole)

1210/DP

Lifting Pole

1292/6/LP

Grab Handles

1210/GRAB

1213/GRAB

1216/GRAB

Backrest

3004

Cromer Cantilever (Bed Cradle)

1369

3-cell Dynamic Mattress

DYN/COM/TRIO

2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

Standard Mattress

MAT/SOFT

Visco Elastic Mattress

MAT/SOFT/VE

Thin Mattress

MAT/SOFT/THIN

Standard Side Rails

UNIVERSAL/SR/MAG

1377/CHROME

High Side Rails

UNIVERSAL/SR/HI/MAG

Side Rail Pads

1316/PS/UNIVERSAL

1316/PS/CREAM

1377/PS

DIMENSIONS:

15cm (H) (Under bed clearance)

100cm (W) x 221cm (L)

BOSTON BED WITH FOUR PLASTIC FEET:

Product Code: 1202/5/PF/PF

BOSTON BED 2 BRAKE AND 2

SWIVEL CASTORS:

Product Code: 1202/5/4/4B

