

BRADSHAW BARIATRIC LOW NURSING CARE BED



The Bradshaw Bariatric Low bed has been designed to provide a robust and sturdy bed for Bariatric use in the Community, Nursing and Residential Care markets. It complements the well established Bradshaw range of beds and provides the opportunity to use a mix of coordinated beds within a caring environment.

It is manufactured with high quality wooden ends and is surprisingly quick and easy to assemble even within limited spaces. The easy assembly is not compromised by the significant increase in engineering specification that has been incorporated into the bed to ensure that it can achieve a maximum user weight of 320kg (50 stone).

The electrically profiling 4 section bed is supplied with a low height, metal slatted mattress platform. There is no provision for use of an integrated lifting pole; instead, for bariatric users we would recommend the use of an appropriate hoist.



KEY FEATURES:

- Low electrically operated variable height for ease of access and to assist the work of the carer, reducing patient lifting and handling
- Low platform height of 26 - 62.5cm
- Maximum user weight of up to 320kg (50 stone)
- Lockable handset controls all electrically operated functions
- Four brake castors, for added safety
- Electrically operated rising backrest and knee-break functions
- Infinitely variable tilt and reverse tilt positioning (Trendelenburg)
- Attractive domestic appearance with high quality wooden head ends, foot ends and full length side rails
- Full range of accessories

ACCESSORIES:

Grab Handle:
1215/GRAB

Bed Extension
(includes longer side rails)
1275/BAR/EX

Cromer Cantilever (Bed Cradle)
1369

2-cell Dynamic Mattress
DYN/BAR

Visco Elastic Mattress
MAT/ACCL/VE/BAR

Visco Elastic Mattress
Extension

MAT/ACCL/VE/BAR/EX

Standard Side Rails

(set of 4)

1275/BAR/SR fitted as standard

Siderail Height Extension (pair)

1275/5/BAR

Side Rail Pads

1316/PS/74/CREAM



SIDE RAIL HEIGHT EXT



BED EXTENSION



MATTRESS EXTENSION



BARIATRIC GRAB HANDLE

DIMENSIONS:

26cm - 62.5cm (Variable platform height) x
136cm (W) x 223.5cm (L)

BRADSHAW BARIATRIC LOW BED:

Product Code: 1275/BAR/LOW