

BRADSHAW LOW NURSING CARE BED



The Bradshaw Low bed has been developed to offer the user the opportunity to set the mattress platform at a very low height of 20cm (8"). It is manufactured with high quality wooden ends, specifically for Community and Nursing Home use, and combines quick and easy assembly, even within limited spaces.

The Bradshaw Low electrically profiling bed, like the Bradshaw standard bed, is particularly versatile and is delivered with wooden side rails and metal slatted mattress platform.



KEY FEATURES:

- Low mattress platform, variable height of 20cm - 60cm
- Electrically operated variable height for ease of access and to assist the work of the carer, reducing patient lifting and handling
- Maximum user weight of up to 180kg (28 stone)
- Electrically operated rising backrest & knee-break functions.
- Infinitely variable tilt and reverse tilt positioning (Trendelenburg)
- Four brake castors, for ease of manoeuvrability and added safety
- Attractive domestic appearance with high quality wooden head ends, foot ends and full length side rails
- Specifically designed, profiling, static & dynamic mattresses available (See Surfaces section Page 40)
- Full range of bed accessories (See Page 36)
- Lockable handset controls all electrically operated functions
- Premium wooden detailing on mattress platform

ACCESSORIES:

Lifting Pole

1275/LP

Grab Handle

1214/GRAB

Bed Extension (includes longer side rails)

1275/EX

Cromer Cantilever (Bed Cradle)

1369

3-cell Dynamic Mattress

DYN/COM/TRIO

2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

Standard Mattress

MAT/SOFT

Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

Visco Elastic Mattress

MAT/SOFT/VE

MAT/ACCL/VE

Thin Mattress

MAT/SOFT/THIN

Standard Mattress Extension

MAT/SOFT/EX

Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

Standard Side Rails (set of 4)

1275/003 fitted as standard

Siderail Height Extension

1275/5

Standard Side Pads

1316/PS/74/CREAM

Siderail Mesh Pads

1275/PS

Crash Mat (Premium)

MAT/FALL

Crash Mat (Basic)

MAT/FALL/BAS

Mesh Cradle

1275/EPS

DIMENSIONS:

20cm - 60cm (Variable platform height)

105.5cm (W) x 222cm (L)

115.5cm (Overall bed height at highest point)

BRADSHAW LOW BED:

Product Code: 1275/M/LOW