



# DURHAM NURSING CARE BED

The Durham bed is specifically designed for Community and Nursing Home use. It has a rising backrest and interchangeable head and foot ends allowing it to blend easily into a domestic environment.



## KEY FEATURES:

- Rising backrest (mechanically assisted)
- Hydraulically operated variable height adjustment (40.5cm - 79cm)
- Maximum user weight of up to 180kg (28 stone)
- Simple to operate foot pedal (optional 2nd pedal for dual side usage)
- Two braked swivel castors, for ease of manoeuvrability and safety
- Specifically designed, static & dynamic mattresses available (See surfaces section - Page 40)
- Full range of bed accessories (See Page 36)

## ACCESSORIES:

### IV (Transfusion Pole)

1210/DP

### Lifting Pole

1292/6/LP

### Grab Handle

1210/GRAB

1213/GRAB

1216/GRAB

### Cromer Cantilever (Bed Cradle)

1369

### 3-cell Dynamic Mattress

DYN/COM/TRIO

### 2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

### Standard Mattress

MAT/SOFT

### Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

### Visco Elastic Mattress

MAT/SOFT/VE

### Thin Mattress

MAT/SOFT/THIN

### Standard Mattress Extension

MAT/SOFT/EX

### Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

### Standard Side Rails

UNIVERSAL/SR/MAG

1377/CHROME

### High Side Rail

UNIVERSAL/SR/HI/MAG

### Side Rail Pads

1316/PS/UNIVERSAL

1377/PS

1316/PS/CREAM

### 2nd Pedal for Durham Bed

1256/PK

## DIMENSIONS:

40.5cm - 79cm (Variable platform height) x  
100cm (W) x 221cm (L)

## DURHAM NURSING CARE BED:

**Product Code: 1258/BED**

