

# INDEPENDENCE INNOV8 RANGE



The new Independence Innov8 range of hospital beds covers every need from a utility ward bed to high dependency requirements.

The beds have been designed from the ground up to satisfy the needs of both the hospital team and the patient. With the input of Tissue Viability Nurses, Infection Control Nurses and Manual Handling specialists, from individual meetings to focus groups, the important features of the beds have been defined, refined and implemented. The key design philosophy, therefore, has been to develop a bed which is easy to use, supports world-class infection control, is reliable and easily maintained, and incorporates superb ergonomics from both a manual handling and patient comfort point of view.

All Independence beds are designed and manufactured in the UK, which brings some very useful benefits. Of course, it means that access to parts, service and maintenance support, and technical expertise is very straightforward, but it goes beyond that.



INNOV8/BE/LC WITH MONITOR SHELF

We are all concerned about our carbon footprint and our impact on the climate, and local manufacture ensures that each bed travels the least distance between the production facility and the hospital, minimising the amount of fuel used to deliver the equipment. Add to this the nightpack option, which uses extremely low-power high-brightness LEDs that can help eliminate the need for more expensive ward lighting, and you have a real ecological winner.

The low bed height helps to reduce the consequences of patient falls, and the wide range of castor options means that you can optimise the balance between bed height, manoeuvrability and tramlining on ridged floor surfaces.


Because of modern infection control requirements the bed has been designed to be exceptionally easy to clean and decontaminate, and to reduce the risk of cross infection. Solid platform panels mean that patient body fluids cannot easily drip onto the bed frame or the floor, and the use of Biocote means that the bed frame itself forms part of the arsenal of infection control processes within the ward or ICU.




INNOV8/2000 WITH OPTIONAL SIDE RAILS

## BED END INFILL COLOUR OPTIONS:

All Innov8 beds feature interchangeable head and foot ends. The standard ends are coloured Aqua, and 5 other colour options are available to help you to manage which beds are associated with which wards.





# INDEPENDENCE INNOV8/2000

The INNOV8/2000 is both the standard ward bed and the base platform to which any accessory can be added. Whatever your specific requirement, you can customise the bed that satisfies your needs by adding the appropriate accessories, giving you total flexibility with regard to bed specification.



For mattresses see pages 40 - 56



## KEY FEATURES:

- Easy to clean and decontaminate
- Incorporates Biocote®, a silver oxide antibacterial coating which helps to reduce the risk of cross infection
- Solid bed panels for maximum infection control and ease of cleaning, with the incorporation of ridges for breathability.
- Electrically operated backrest, height and knee-break
- Infinitely variable electrically operated Trendelenburg and reverse Trendelenburg positions
- Auto-contour - simultaneous adjustment of backrest and knee-break
- Auto-regression, avoids surface pinching or patient sliding when in Fowler position
- 193mm underbed clearance for hoist access
- Fitted with 12.5cm single castors
- Safe working load of up to 255kg (40 stone)
- Very low bed height to minimise consequences of patient falls
- Battery backup
- Patient Handset and Nurse Control Panel, with handset lockout (via the Nurse Control Panel)
- Dual sided CPR levers
- IP54 standard ingress protection
- Very easy to manoeuvre and work with, reducing back or handling problems
- Removable head/foot ends with colour coding options for estates management
- UK manufacture and support
- Modular, hence easy to maintain
- Total patient comfort
- Simple to use

## FACTORY FITTED ACCESSORIES:

Night pack INNOV8/NP	15cm single castors INNOV8/SC150
Standard cantilever side rails INNOV8/SR	20cm single castors INNOV8/SC200
High cantilever side rails INNOV8/HSR	12.5cm double castors INNOV8/DC125
X-ray backrest INNOV8/XB	15cm double castors INNOV8/DC150
IP 66 INNOV8/IP66	

## GENERAL ACCESSORIES:

Lifting pole INNOV8/LP	Monitor shelf INNOV8/MS
Drip pole INNOV8/DP	Oxygen Cylinder Carrier INNOV8/OC
Traction – Balkan Beam INNOV8/TBB	
Traction – Foot End INNOV8/TFE	


Note: Minimum quantities apply for IP66, special castors and bed end infill I colours.


## DIMENSIONS:

40.7cm - 80.7cm (H) x 101.6cm (W) x 213cm (L)  
(233cm fully extended)

## INDEPENDENCE INNOV8/2000:

Product Code: **INNOV8/2000**





## INDEPENDENCE INNOV8/4000

The INNOV8/4000 is an advanced ward bed, featuring split side rails, simplicity and ease of use, excellent infection control, and superb ergonomics. As with all INNOV8 beds the wide range of options allows for almost unlimited customisation to suit individual hospital needs.



For mattresses see pages 40 - 56



### KEY FEATURES:

- Split side rails offering easy and safe access to patients
- Easy to clean and decontaminate
- Incorporates Biocote®, a silver oxide antibacterial coating which helps to reduce the risk of cross infection
- Solid bed panels for maximum infection control and ease of cleaning, with the incorporation of ridges for breathability.
- Electrically operated backrest, height and knee-break
- Infinitely variable electrically operated Trendelenburg and reverse Trendelenburg
- Auto-contour - simultaneous adjustment of backrest and knee-break
- Auto-regression, avoids surface pinching or patient sliding when in Fowler position
- 193mm underbed clearance for hoist access
- Fitted with 12.5cm single castors
- Safe working load of up to 255kg (40 stone)
- Very low bed height to minimise consequences of patient falls
- Battery backup
- Patient Handset and Nurse Control Panel, with handset lockout (via the Nurse Control Panel)
- Dual sided CPR levers
- IP54 standard ingress protection
- Very easy to manoeuvre and work with, reducing back or handling problems
- Removable head/foot ends with colour coding options for estates management
- UK manufacture and support
- Modular, hence easy to maintain
- Total patient comfort
- Simple to use

### FACTORY FITTED ACCESSORIES:

Night pack  
INNOV8/NP  
X-ray backrest  
INNOV8/XB  
IP 66  
INNOV8/IP66  
15cm single castors  
INNOV8/SC150

20cm single castors  
INNOV8/SC200  
12.5cm double castors  
INNOV8/DC125  
15cm double castors  
INNOV8/DC150

Note: Minimum quantities apply for IP66, special castors and bed end infill I colours.

### GENERAL ACCESSORIES:


Lifting pole  
INNOV8/LP  
Drip pole  
INNOV8/DP  
Traction – Balkan Beam  
INNOV8/TBB  
Traction – Foot End  
INNOV8/TFE  
Oxygen cylinder carrier  
INNOV8/OC  
Monitor shelf  
INNOV8/MS

### DIMENSIONS:

40.7cm - 80.7cm (H) x 101.6cm (W) x 213cm (L)  
(233cm fully extended)

### INDEPENDENCE INNOV8/4000:

Product Code: **INNOV8/4000**



## INDEPENDENCE INNOV8/6000

The advanced INNOV8/6000 bed is designed for high dependency patients. It incorporates split side rails, an X-ray translucent backrest and an oxygen cylinder carrier. Key accessories include a monitor shelf.



The bed has been designed with wound care management as a major priority, yet is still exceptionally easy to use. The X-ray backrest reduces the need for patient transfers to alternative platforms, improving patient outcomes and reducing manual handling of patients who are often not mobile or able themselves to provide any load bearing support or assistance.



For mattresses see pages 40 - 56



### KEY FEATURES:

- Split side rails offering easy and safe access to patients
- X-ray translucent backrest
- O<sub>2</sub> cylinder carrier
- Easy to clean and decontaminate
- Incorporates Biocote®, a silver oxide antibacterial coating which helps to reduce the risk of cross infection
- Solid bed panels for maximum infection control and ease of cleaning, with the incorporation of ridges for breathability.
- Electrically operated backrest, height and knee-break
- Infinitely variable electrically operated Trendelenburg and reverse Trendelenburg positions
- Auto-contour - simultaneous adjustment of backrest and knee-break
- Auto-regression, avoids surface pinching or patient sliding when in Fowler position
- Fitted with 12.5cm single castors
- Safe working load of up to 255kg (40 stone)
- 193mm underbed clearance for hoist access
- Very low bed height to minimise consequences of patient falls
- Battery backup
- Patient Handset and Nurse Control Panel, with handset lockout (via the Nurse Control Panel)
- Dual sided CPR levers
- IP54 standard ingress protection
- Very easy to manoeuvre and work with, reducing back or handling problems
- Removable head/foot ends with colour coding options for estates management
- UK manufacture and support
- Modular, hence easy to maintain
- Total patient comfort
- Simple to use

### FACTORY FITTED ACCESSORIES:

Night pack INNOV8/NP	20cm single castors INNOV8/SC200
IP 66 INNOV8/IP66	12.5cm double castors INNOV8/DC125
15cm single castors INNOV8/SC150	15cm double castors INNOV8/DC150

### GENERAL ACCESSORIES:

Lifting pole INNOV8/LP	Traction – Foot End INNOV8/TFE
Drip pole INNOV8/DP	Monitor shelf INNOV8/MS
Traction – Balkan Beam INNOV8/TBB	

Note: Minimum quantities apply for IP66, special castors and bed end infill I colours.

### DIMENSIONS:

40.7cm - 80.7cm (H) x 101.6cm (W) x 213cm (L)  
(233cm fully extended)

### INDEPENDENCE INNOV8/6000:

Product Code: **INNOV8/6000**

Tel: +44 (0)1422 233000 sales@sidhil.com

||

www.sidhil.com Fax: +44 (0)1422 233010