



INDEPENDENCE VERSA



The Independence Versa is a fully profiling, electrically operated hospital bed, developed specifically to meet the varying needs of the acute sector. It achieves this by being versatile, comfortable and user friendly for both the patient and the carer.

With its ergonomically designed longer backrest, active seat system and 'easy-clean' design it has responded to the ever-changing needs of today's healthcare environment.



KEY FEATURES:

- Bed lengthening facility to cater for taller patients
- 'Active-Seat' to reduce 'shearing' and smooth the transfer to seat position
- Standard and High side rail options, fully compliant with current acute nursing requirements
- 'Easy-Use' locking side rail mechanism
- Integrated X-ray backrest panel (factory fitted option)
- Handset lockout box for easy bed control management
- Maximum user weight of up to 267kg (42 stone)
- Ventilated, removable, easy to clean mattress platform sections
- Range of pressure relieving, foam and alternating dynamic therapy mattress systems developed for use with the Independence bed.
- 3 year warranty (1st year full parts & labour, 2 further years parts only)

ACCESSORIES:

Oxygen Cylinder Carrier
2250/OX

Lifting Pole, Strap & Handle -
1292/6/LP/W

X-Ray Backrest Panel -
(specify at time of order)
2260/AQ/X

Range of Static & Dynamic
Mattresses - specifically designed for
use with Independence beds
(see surfaces section)

High or Standard Side Rails -
2200/RSR/HI/CH (High)
2200/RSR/2/CH (Standard)
IV (Transfusion) Pole -
2200/DP

DIMENSIONS:

38cm - 76.5cm (Mattress platform height) x
100.7cm (W) x 212cm - 232cm (L)
190.6cm - 216cm (Internal length)

INDEPENDENCE VERSA:

Product Code: 2260/AQ (Side rails not included)

