


# INSPIRATION COT



Designed in conjunction with leading UK Paediatric specialists, the Inspiration children's cot is colourful and aesthetically pleasing whilst still featuring innovative components. To maximise patient care, up to two removable transfusion poles can be fitted to the cot. Additionally, a dual sided, linked, tracking and braking system has been added for safety.



LULLABY VE



HINGED ACCESS

## KEY FEATURES:

- Electrically operated
- Variable height
- Two-way tilting mattress platform
- Uniquely designed drop-down telescopic side rails
- Optional manual CPR mechanism
- Maximum user weight when not in tilt: 160kg (25 Stone)
- Maximum user weight when in tilt: 50kg (8 Stone)
- Removable head and foot ends give flexible patient access (combined with the drop-down telescopic side rails allows four sided access)
- Specially designed visco-elastic mattress included as standard

## ACCESSORIES:

Lullaby VE  
MAT/LULL/VE  
Lullaby Pressure Cradle  
DYN/LULL  
Cot Traction System  
TRACT/PAED

Extra Transfusion Poles  
9000/DP  
Oxygen Cylinder Carrier  
9000/OX  
Lockout Box  
9000/LOCKOUT



COT TRACTION SYSTEM

## DIMENSIONS:

58.5cm - 92cm (Mattress platform height)  
82cm (W) x 145cm (L)

## INSPIRATION COT INCLUDES

### LULLABY VE MATTRESS:

**Product Code: 9000**

**9000 COT WITH HINGED ACCESS AND CPR SYSTEM INCLUDES LULLABY VE MATTRESS:**

**Product Code: 9003**