



SOLITE LASER 2 BED



The Solite Laser 2 bed has been developed to offer all the advantages you can expect from the Vari-Height bed (see page 19), together with an electrically operated backrest. The backrest offers further handling benefits to the carer and a higher degree of independence to the user.

The bed has been designed for quick and simple, yet safe assembly by one person working alone, without the need for special hand tools. Adding to ease of storage, movement and assembly, the bed is made up of four component sections, the heaviest being only 22kg.



The 12.7cm Solite Bed Extension can be used in conjunction with a mattress extension

On the Solite Laser range of beds, the mattress platform is capable of lowering to a height of only 31.5cm and can be raised to a height of 72cm using a personal handset (from floor to top of mattress platform).

The beds are all capable of taking a maximum user weight of 180kg (28 stone).

ACCESSORIES:

IV (Transfusion Pole)

1210/DP

Lifting Pole

SOLITE/LASER/LP

Grab Handle

1210/GRAB

1213/GRAB

1216/GRAB

Bed Extension

SOLITE/LASER/EXT

Cromer Cantilever (Bed Cradle)

1369

3-cell Dynamic Mattress

DYN/COM/TRIO

2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

Standard Mattress

MAT/SOFT

Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

Visco Elastic Mattress

MAT/SOFT/VE

MAT/ACCL/VE

Thin Mattress

MAT/SOFT/THIN

Standard Mattress Extension

MAT/SOFT/EX

Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

Mattress Infills

SOLITE/SR/INFILL

MAT/INFILLS

Standard Side Rails

SOLITE/SR/2

UNIVERSAL/SR/MAG

1377/CHROME

High Side Rail

UNIVERSAL/SR/HI/MAG

Side Rail Pads

1316/PS/CREAM

1377/PS

1316/PS/UNIVERSAL

DIMENSIONS:

31.5cm - 72cm (Variable platform height) x

94.5cm (W) x 219cm (L)

SOLITE LASER 2 BED WITH RISING BACKREST:

Product Code: SOLITE/LASER/2

