



SOLITE LASER 4 BED

The Solite Laser 4 electrically operated, fully profiling bed has been developed specifically for use in the home care and nursing care sectors. It delivers a high degree of independence and freedom of movement through profiling features for users and carers alike. A lockout handset is now included as standard. This allows the carer to limit the use of any or all of the functions on the beds handset.

It has been designed for quick and simple, yet safe, assembly by one person working alone, adding to ease of storage, movement and assembly. The Solite Laser 4 bed has a battery backup option which enables a limited amount of bed movements in the case of mains power supply failure.



KEY FEATURES:

- Electrically operated variable height for ease of access and to assist the work of the carer, reducing patient lifting and handling
- Maximum user weight of up to 180kg (28 stone)
- Single handset controls all electrically operated functions
- Electrically operated rising backrest section and knee break section
- Infinitely variable tilt and reverse tilt positioning (Trendelenburg)
- Four brake castors, for ease of manoeuvrability and added safety
- Attractive domestic appearance with light oak head and foot ends, and special scratch-resistant magnolia coating to all metal parts
- Specifically designed profiling, static & dynamic mattresses available (See Surfaces section - Page 40)
- Full range of bed accessories (See Page 36)

ACCESSORIES:

IV (Transfusion Pole)
I210/DP

Lifting Pole

SOLITE/LASER/LP

Grab Handle

I210/GRAB

I213/GRAB

I216/GRAB

Bed Extension

SOLITE/LASER/EXT

Cromer Cantilever (Bed Cradle)

I369

3-cell Dynamic Mattress

DYN/COM/TRIO

2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

Standard Mattress
MAT/SOFT

Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

Visco Elastic Mattress

MAT/SOFT/VE

MAT/ACCL/VE

Thin Mattress

MAT/SOFT/THIN

Solite Laser Bed extension

SOLITE/LASER/EXT

Standard Mattress Extension

MAT/SOFT/EX

Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

Mattress Infills

UNIVERSAL/SR/WH

Standard Side Rails

SOLITE/SR/2

UNIVERSAL/SR/MAG

I377/CHROME

High Side Rails

UNIVERSAL/SR/HI/MAG

Side Rail Pads

I316/PS/CREAM

I377/PS

I316/PS/UNIVERSAL

Spare Lockable Handset

CONTROL/1038

Handset Locking Key

CONTROL/049

DIMENSIONS:

31.5cm - 72cm (Variable platform height) x
94.5cm (W) x 219cm (L)

SOLITE LASER 4 FULLY PROFILING BED:

Product Code: SOLITE/LASER/4

SOLITE LASER 4 FULLY PROFILING BED

WITH BATTERY BACKUP:

Product Code: SOLITE/LASER/4/BB

