



SOLITE LASER VARI-HEIGHT BED



The Solite Laser Vari-Height has been developed to replace the traditional, manually adjustable, variable height bed. Through the addition of electrical actuators to raise and lower the bed and utilising a domestic design approach, the Solite Laser Vari-Height is an ideal, general purpose, easily moveable, community care bed.



On the Solite Laser range of beds, the mattress platform is capable of lowering to a height of only 31.5cm and can be raised to a height of 72cm using a personal handset (from floor to top of mattress platform).

The beds are all capable of taking a maximum user weight of 180kg (28 stone).

• Please note no lifting pole is available for the Solite Laser Vari-Height

ACCESSORIES:

IV (Transfusion Pole)

1210/DP

Grab Handle

1210/GRAB

1213/GRAB

1216/GRAB

Bed Extension

SOLITE/LASER/EXT

Backrest

3004

Cromer Cantilever (Bed Cradle)

1369

3-cell Dynamic Mattress

DYN/COM/TRIO

2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

Standard Mattress

MAT/SOFT

Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

Visco Elastic Mattress

MAT/SOFT/VE

MAT/ACCL/VE

Thin Mattress

MAT/SOFT/THIN

Standard Mattress Extension

MAT/SOFT/EX

Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

Mattress Infills

SOLITE/SR/INFILL

MAT/INFILLS

Standard Side Rails

SOLITE/SR/2

UNIVERSAL/SR/WH

I377/CHROME

High Side Rail

UNIVERSAL/SR/HI/WH

Side Rail Pads

I316/PS/CREAM

I377/PS

I316/PS/UNIVERSAL

DIMENSIONS:

31.5cm - 72cm (Variable platform height) x
94.5cm (W) x 219cm (L)

SOLITE LASER VARI-HEIGHT BED:

Product Code: SOLITE/LASER/VH

