



## YORK NURSING CARE BED

The York bed is specifically designed for community and nursing home use. It has interchangeable head and foot ends allowing it to blend easily into a domestic environment.



### KEY FEATURES:

- Hydraulically operated variable height adjustment (40.5cm - 79cm)
- Maximum user weight of up to 180kg (28 stone)
- Simple to operate foot pedal (optional 2nd pedal for dual side usage)
- Two braked swivel castors, for ease of manoeuvrability and safety
- Specifically designed, static & dynamic mattresses available (see surfaces section - Page 40)
- Full range of bed accessories (See Page 36)

### ACCESSORIES:

#### IV (Transfusion Pole)

1210/DP

#### Lifting Pole

1292/6/LP

#### Grab Handles

1210/GRAB

1213/GRAB

1216/GRAB

#### Backrest

3004

#### Cromer Cantilever (Bed Cradle)

1369

#### 3-cell Dynamic Mattress

DYN/COM/TRIO

#### 2-cell Dynamic Mattress

DYN/COM/PLU/2

DYN/COM/SOLO/2

#### Standard Mattress

MAT/SOFT

#### Profiling Mattress

MAT/SOFT/CON

MAT/ACCL/PRO

#### Visco Elastic Mattress

MAT/SOFT/VE

#### Thin Mattress

MAT/SOFT/THIN

#### Standard Mattress Extension

MAT/SOFT/EX

#### Visco Elastic Mattress Extension

MAT/ACCL/VE/EX

#### Standard Side Rails

UNIVERSAL/SR/MAG

1377/CHROME

#### High Side Rail

UNIVERSAL/SR/HI/MAG

#### Side Rail Pads

1316/PS/UNIVERSAL

1316/PS/CREAM

1377/PS

#### 2nd Pedal for York Bed

1256/PK

### DIMENSIONS:

40.5cm - 79cm (Variable platform height)  
x 100cm (W) x 221cm (L)

### YORK NURSING CARE BED:

Product Code: 1256/BED

