



CAMBOURNE VARIABLE HEIGHT COUCH



KEY FEATURES:

- Suitable for a wide range of uses
- Hydraulic variable height couch
- 4 brake/swivel castors, for ease of manoeuvrability and safety
- Maximum User Weight 180kg (28 stone)
- Durable, white epoxy coating with Biocote®

ACCESSORIES:

2nd Pedal for Cambourne Couch (not pictured)
3010/PK

'Swing-Away' side-rails (not pictured)
1348/SR

OVERALL DIMENSIONS:

46cm - 83cm (H) x 60cm (W) x 186cm (L)
Product Weight: 70kg

CAMBOURNE VARIABLE HEIGHT COUCH:

Product Code: 3010BLK (Black)
Product Code: 3010TAN (Tan)

SINGLE & DOUBLE COUCH STEPS

These attractively designed steps, which are available in either silver or white paint finish with black treads, have been developed to retain the weight loading and grooved anti-slip rubber treads previously incorporated in Doherty couch steps.

Additionally, a new feature has been added to the products through the introduction of adjustable height points within the feet to aid stability on uneven surfaces. The base of the feet also includes an anti-slip design.

DOUBLE TIER STEP :

Product Code: SIGSTEP2/SILVER
Product Code: COU78/WHITE

SINGLE TIER STEP:

Product Code: SIGSTEP1/SILVER
Product Code: COU73/WHITE



User weight: 160kg (25 stone)

