


MARLETT FIXED HEIGHT COUCH



Marlett Couch with multi-drawer & cupboard.
Finished in maple veneer with soft bevelled
edges and silver epoxy coating.



**MARLETT COUCH WITH
MULTI-DRAWER & CUPBOARD:**
Product Code: **SIGCOU10/(Colour)**

KEY FEATURES:

- 4 drawer storage unit and cupboard with adjustable shelf
- Integral pull-out work surface in both cupboard and drawer unit
- Multi-position adjustable backrest
- Maximum user weight 160kg (25 stone)
- Fire retardant upholstery
- Available in 12 modern colour options

ACCESSORIES:

Lithotomy Stirrups & Leg Holders (not pictured)
SIGCOU14

Replacement Holders for Stirrups (not pictured)
SIGCOU4



MARLETT COUCH WITH MULTI-DRAWER & SINGLE LEG:
As above but with four drawer storage unit
Product Code: **SIGCOU11/(Colour)**



MARLETT COUCH WITH CUPBOARD & SINGLE LEG:
As above but with cupboard unit with adjustable shelf
Product Code: **SIGCOU12/(Colour)**



OVERALL DIMENSIONS:
83cm (H) x 60cm (W) x 183cm (L)

