


PURLEY FIXED HEIGHT COUCH



Purley Couch (with optional drawer units). Strong and stable couch frame construction with a fixed seat height of 76cm (30"). Available in 12 modern colour options.



KEY FEATURES:

- Self-locking, multi-position adjustable backrest
- Maximum user weight of up to 160kg (25 stone)
- 5cm thick foam top for patient comfort
- Paper roll holder (paper roll not supplied)
- Durable, white epoxy coating with Biocote®
- Available in 12 modern colour options

ACCESSORIES:

- Drawer Unit (per drawer)
COU96
- Lithotomy Stirrups & Leg Holders
(not pictured)
COU71
- Replacement Holders for Stirrups
(not pictured)
COU76



OVERALL DIMENSIONS:
76cm (H) x 60cm (W) x 183cm (L)
Product Weight: 34kg

PURLEY COUCH:
Product Code: **COU31/(Colour)**

