



SALISBURY

ELECTRICALLY CONTROLLED, MULTI-FUNCTION VARIABLE HEIGHT COUCH



The Salisbury, electrically controlled, multi-function, variable height couch has been designed to meet the ever growing needs of the busy practice. This design has evolved from our wide ranging experience serving the needs of plinth, couch and bed customers, in addition to extensive user group feedback.



KEY FEATURES:

- Electric hand-set control
- 'Swing-Away' side-rails fitted as standard
- Maximum user weight 160kg (25 stone)
- 4 brake/swivel castors, for ease of manoeuvrability and safety
- Adjustable head/back-rest section with optional breathing hole
- Comfortable foam surface
- Available in 12 modern colour options
- 3 year warranty (1st year parts and labour, 2 further years parts only)

ACCESSORIES:

- Complete Lithotomy Set (including brackets, removable straps and poles) COU34
- Couch End (dual use - push/steer bar & patient safety feature) COU33
- IV/Transfusion Pole (not pictured) COU35



LITHOTOMY SET



COUCH END



OVERALL DIMENSIONS:

52cm - 85cm (H) x 61cm (W) x 185cm (L)
Product Weight: 70kg

SALISBURY COUCH:

Product Code: COU32 (Colour)

SALISBURY COUCH WITH BREATHING HOLE (not pictured):

Product Code: COU32/BH/(Colour)

