



# VERDANA

## VARIABLE HEIGHT COUCH

Four section electrically operated variable height treatment couch.



### KEY FEATURES:

- Split leg section
- Breathing hole (optional)
- Height range 46.5cm - 93cm
- Maximum user weight 160kg (25 stone)
- Manually adjustable self-locking backrest and foot section
- Head down tilt and feet-up positions
- Electric hand set control
- Retractable castor for ease of movement
- Simple levelling for uneven floors
- Attractive silver epoxy coating with Biocote®
- Available in 12 modern colour options
- Paper roll holder
- 3 year warranty (1st year parts and labour, 2 further years parts only)

### ACCESSORIES:

Arm Rests (pair)  
SIGCOU18/(COLOUR)

Replacement Holders for Stirrups (pair)  
SIGCOU4

Side Rails (pair)  
SIGCOU19

Lithotomy Stirrups & Leg Holders (pair)  
SIGCOU21



### OVERALL DIMENSIONS:

46.5cm - 93cm (H) x 65cm (W) x 189cm (L)  
Product Weight: 83kg

### VERDANA COUCH:

Product Code: **SIGCOU27/(Colour)**

### VERDANA COUCH:

As above with breathing hole

Product Code: **SIGCOU27/BH/(Colour)**

