



# WINCHESTER COUCH DRAWER SYSTEM



The Doherty Winchester Couch Drawer System has been designed to complement the well established Winchester Couch. The system allows the user to select numerous combinations of the three core items.



### KEY FEATURES:

- From only three separate products an almost infinite range of combinations is possible
- The system comprises a drawer, shelf and trolley base which can be combined to suit individual requirements
- All three components are available to purchase separately and are easily fitted to the couch or trolley base as required
- Attractive and durable mild steel epoxy coated finish with Biocote®

**SINGLE DRAWER** - Product Code: **COU92/DRW**

**SHELF** - Product Code: **COU92/SHF**

**TROLLEY BASE** - Product Code: **COU92/TLY**  
(Unit fits onto bottom drawer to give wheeled mobility. Does not include drawers)



3 single drawers arranged vertically and one shelf unit



7 single drawers, 1 shelf unit and 1 trolley base: 3 arranged vertically with shelf unit, and 4 attached to trolley base unit



3 single drawers and one shelf unit arranged horizontally

	WEIGHT	DEPTH	WIDTH	HEIGHT
DRAWER	11.5kg	66cm	56.5cm	13cm
SHELF	8kg	66cm	56.5cm	8.5cm
TROLLEY BASE	6kg	65cm	56.5cm	11cm

