


WINCHESTER FIXED HEIGHT COUCH



The Doherty Winchester Couch has been designed to provide an exceptionally robust, sturdy couch with a modern look and feel at a competitive price. The needs of the modern practice have been fully considered from the General Practitioner and Nursing Staff through to patient comfort.

The couch is designed to coordinate with the Doherty Howarth trolley range and the Doherty solid panel screen system.



KEY FEATURES:

- Unique paper roll holder supplied as standard that can be fitted either within the couch frame or externally
- Easy to clean, hygienic, epoxy coated white frame
- Maximum user weight 160kg (25 stone)
- Extremely reliable self-locking multi-position adjustable backrest
- 50mm thick foam top for patient comfort
- Unstitched, rounded upholstery corners for hygiene and cleanliness
- Available in 12 modern upholstery colour options
- Simple levelling for uneven floors
- Sculptured tube frame to ease back rest lifting


ACCESSORIES:

Complete Lithotomy Set including length adjustable brackets, poles and removable straps
COU91

Replacement Holders for Stirrups
COU76



Self locking multi-position adjustable backrest



Adjustable lithotomy bracket



Paper roll holder



Option of internal or external fitting for paper roll holder



OVERALL DIMENSIONS:

80cm (H) x 70.5cm (W) x 184cm (L)
paper roll holder "inboard"
80cm (H) x 70.5cm (W) x 188cm (L)
paper roll holder "outboard"
Product weight (without drawers): 40kg

WINCHESTER FIXED HEIGHT COUCH

Product Code: **COU90/(colour)**

