



DOHERTY 3 SECTION PLINTH

ELECTRICALLY CONTROLLED, VARIABLE HEIGHT



The Doherty 2 and 3 Section Plinths have been developed specifically for applications throughout the modern GP practice or Physiotherapy environment. The Plinths have been designed to a very high specification for optimum flexibility. Available in 12 upholstery colour options, they are designed to coordinate with other products within the extensive Doherty product range, including the Howarth trolleys.



KEY FEATURES:

- Electrically operated variable height range from 46cm to 85cm
- Easy lift gas assisted backrest and foot section
- Capable of both head down tilt and feet up positions
- Unique self levelling system when lowering plinth
- Electric handset control for height adjustment
- Maximum user weight 160kg (25 stone)
- Breathing hole and paper roll holders fitted as standard
- Distinctive Doherty style upholstery with rounded corners for hygiene and cleanliness
- Twin retractable castors for ease of movement
- Available in 12 modern colour options
- Attractive white epoxy coating with Biocote®
- Delivered fully assembled
- 3 year warranty for peace of mind (1st year parts and labour, 2 further years parts only)
- Comprehensive range of stylish accessories available
- Handset holder
- Dual sided operation
- Simple levelling for uneven floors

For accessories please see page 64



OVERALL DIMENSIONS:

46cm - 85cm (H) x 60cm (W) x 193cm (L)
Product Weight: 78kg

DOHERTY 3 SECTION PLINTH:

Product Code: PLE02/(Colour)

Tel: +44 (0)1422 233000 sales@sidhil.com

63

www.sidhil.com Fax: +44 (0)1422 233010

