

# FIRST AID



### SIDHIL TRANSIT CHAIR RANGE

Ideal for use in First Aid rooms, ambulances, retail, industrial and general workplace environments

- Green, wipe-clean PVC chair cover
- Safety belt fitted as standard
- Wide range of options to cover all needs

### DIMENSIONS:

92.5cm (H) x 46cm (W)  
 64.5cm (Length 2 wheel), with footrest 82.5cm  
 66cm (Length 4 wheel), with footrest 82.5cm  
 Maximum user weight: 115kg (18st)

#### 1. TRANSIT CHAIR (2 rear wheels)

Product Code: TRA01

#### 2. TRANSIT CHAIR (2 rear wheels & footrest)

Product Code: TRA02

#### 3. TRANSIT CHAIR (4 wheels)

Product Code: TRA03

#### 4. TRANSIT CHAIR (4 wheels & footrest, front braked)

Product Code: TRA07

#### TRANSIT CHAIR (4 wheels & footrest)

Product Code: TRA04 (main picture)

#### SPARE SAFETY BELT FOR TRANSIT CHAIR

(not pictured)

Product Code: TRA05

#### SPARE PVC SEAT COVER FOR TRANSIT CHAIR

Product Code: TRA06



### CHARD PORTABLE STRETCHER

- Lightweight design - only 6kg
- Folds flat for easy storage and transportation
- Aluminium poles with rubber hand grips
- Aluminium feet prevent the stretcher bed touching floor
- Optional carry case (STR57)

### DIMENSIONS:

14.5cm (H) x 58.5cm (W) x Open: 216.5cm (L)  
 Folded: 9cm (H) x 17cm (W) x 107.5cm (L)  
 Maximum user weight: 115kg (18st)

Product Code: STR01

### CARRY CASE FOR CHARD PORTABLE

### STRETCHER DIMENSIONS:

30cm (W) x 113cm (L)  
 Product weight 0.7kg

Product Code: STR57

