

## SUPPORT PRODUCTS



### PERCHING STOOL (with arms & backrest)

- Designed with seat front 50mm lower than seat back
- Maximum user weight 160kg (25 stone)
- Splayed legs for extra stability
- Extending legs for height adjustment
- Mounted on non-slip rubber feet
- Fixed arms and back support
- Stackable

**Product Code: 1398/EX/WHITE**

### PERCHING STOOL (with arms & backrest)

(As above with non-padded arms & padded backrest)

**Product Code: 1399/EX/WHITE**

### DIMENSIONS:

Seat Front Height 60cm

Seat Area 34.6cm x 29.6cm

Arm Height Above Seat 13.5cm

Width Between Arms 46cm

Leg Height Range 15cm



### WHITBY LIFTING POLE

- Combined lifting pole and standing aid
- Maximum patient upper body user weight - 80kg (13 stone)
- Adjustable strap and handle

### DIMENSIONS:

190cm (H) x 77cm (W)

71cm (Length of arm reach)

102cm (Support leg length)

**Product Code: 1385/1**

### STRAP & HANDLE

Spare for Whitby lifting pole (1385/1)

**Product Code: 1209/032/033**

### ROMFORD PADDED LEG REST

(Must not be used as a step stool)

- Padded top cover as standard with flame retardant vinyl
- Mounted on non-slip rubber feet
- Adjustable height option on each leg
- Height 35.5cm - 50cm

**Product Code: 1375**

### ROMSEY PADDED LEG REST (with four castors)

(Must not be used as a step stool)

- Padded top cover as standard with flame retardant vinyl
- Mounted on four free running castors
- Adjustable height option on each leg
- Height 38.2cm - 52.7cm

**Product Code: 1375/C**

