

## OVERBED TABLES, OVERCHAIR TABLES & CABINETS



The Sidhil range of Overbed Tables have been designed for hospital, nursing, residential care home and community use. All 'Easi-Riser' products are available in both standard or two-piece tilting tops, with a choice of bases to suit different environments and patient requirements.



**EASI-RISER OVERBED TABLE (standard base)**



**DIMENSIONS:**  
72cm - 110cm (H) x 93cm (W) x 40cm (D)  
**Product Code: 3012**

**EASI-RISER OVERBED TABLE (wheel chair base)**



As 3012 with a curved base

**DIMENSIONS:**  
72cm - 110cm (H) x 93cm (W) x 40cm (D)  
**Product Code: 3013**

**KEY FEATURES:**

- Height adjustable from 72cm to 110cm
- Height adjusting handle
- Anti-trapping safety feature
- Fire retardant, light oak laminate top
- Raised upstand
- Four swivel castors
- Magnolia epoxy coating with Biocote®
- Supplied fully assembled

**EASI-RISER OVERBED TABLE**  
(tilting top / standard base)



As 3012 with a tilting top  
• Tilting top adjusts from 0° to 50°


**DIMENSIONS:**  
72cm - 110cm (H) x 93cm (W) x 40cm (D)  
**Product Code: 3014**


**EASI-RISER OVERBED TABLE**  
(tilting top / wheel chair base)



As 3012 with a tilting top and curved base  
• Tilting top adjusts from 0° to 50°

**DIMENSIONS:**  
72cm - 110cm (H) x 93cm (W) x 40cm (D)  
**Product Code: 3015**





## OVERBED TABLES, OVERCHAIR TABLES & CABINETS



**CENTENARY OVERBED TABLE** (self assembly)

- Height adjustable from 62cm to 94cm
- Class 3 fire retardant, light oak laminate top
- Metal frame base
- Tilting top adjusts from 0° to 70°
- Magnolia epoxy coating
- Raised upstands

**DIMENSIONS**

62cm - 94cm (H) x 60cm (W) x 40cm (D)

**Product Code: I337/MAG**

**THAMESMEAD OVERBED TABLE** (with four castors)

As I337/MAG with four castors (easily assembled to base)



**DIMENSIONS** (from floor):

68cm - 100cm (H) x 60cm (W) x 40cm (D)

**Product Code: I337/C/MAG**

**BEDSIDE CABINET**



- Saloon doors for easier access
- Integral towel rail holder at rear
- Four swivel castors
- Fire retardant light oak laminate

**DIMENSIONS**

95.5cm (H, including castors) x 45cm (W) x 45cm (D)

**Product Code: 3021**

**ASHBOURNE OVERCHAIR TABLE**





- Extra-wide overchair table
- Manual height adjustment from 59cm to 83cm
- Raised upstands
- Four swivel castors
- Fire retardant, laminate top
- Magnolia epoxy coating with Biocote®
- Supplied fully assembled

**DIMENSIONS**

59cm - 83cm (H) x 84cm (W) x 36cm (D)

**Product Code: OCT02**





## DRUG CABINETS



The Doherty range of Controlled Drug Cabinets is manufactured in 14 gauge mild steel with a white, hard-wearing epoxy finish. They have been designed to meet the requirements of the Misuse of Drugs (Safe Custody) Regulations 1973 and are ideally suited to the Industrial First Aid and GP Surgery markets.

The design incorporates a composite frame within the cabinet door, together with unique heavy duty internal lock cover-plates, preventing removal of the lock. Heavy duty flanges around the front opening also give excellent security protection.

### WALL UNIT, 1 SHELF, 1 DOOR & 1 LOCK

25cm (H) x 30cm (W) x 22.5cm (D)

**Product Code: CDC1000W**



### WALL UNIT, 1 SHELF, 1 DOOR & 1 LOCK

45cm (H) x 35cm (W) x 15cm (D)

**Product Code: CDC1005W**



### WALL UNIT, 1 SHELF, 1 DOOR & 1 LOCK

40cm (H) x 30cm (W) x 25cm (D)

**Product Code: CDC1009W**



### WALL UNIT, 2 SHELVES, 1 DOOR & 1 LOCK

60cm (H) x 60cm (W) x 25cm (D)

**Product Code: CDC1010W**



### WALL UNIT, 3 SHELVES, 1 DOOR & 2 LOCKS

90cm (H) x 60cm (W) x 25cm (D)

**Product Code: CDC1015W**



### FLOOR UNIT, 3 SHELVES, 1 DOOR & 2 LOCKS

80cm (H) x 60cm (W) x 45cm (D)

**Product Code: CDC1020F**

