


DOHERTY PHLEBOTOMY AND TREATMENT CHAIRS



The new Doherty Phlebotomy and Treatment Chairs have been developed specifically for many applications throughout the modern GP practice, Treatment Room and Hospital environment. Our highly experienced engineering team have developed the specification and executed the design following extensive consultation with specialists around the UK. The range comprises 3 products; a Phlebotomy Chair, a Variable Height Treatment Chair and a Fixed Height Treatment Chair. All 3 products, which are available in 12 modern upholstery colour options, are designed to coordinate with other products (including the Howarth trolleys) within the extensive Doherty product range.

DOHERTY PHLEBOTOMY CHAIR



KEY FEATURES:

- Electrically operated variable height range
- Electric handset control for height adjustment
- Breathing hole available as an optional extra
- The Phlebotomy Chair has the added advantage of also being suitable for use as a plinth. Where cost or space constraints limit the amount of equipment that can be purchased this facility makes the chair particularly attractive
- A paper roll holder is available as an "easi fit" accessory
- Easy lift gas assisted backrest and foot section
- Maximum user weight 160kgs (25 stone)
- Dual sided operation
- Distinctive Doherty style upholstery with rounded corners for hygiene and cleanliness
- Twin retractable castors for ease of movement
- Simple levelling mechanism for uneven floors
- Attractive white epoxy coating with Biocote
- Delivered fully assembled
- Available in 12 modern upholstery colour options
- 3 year warranty for peace of mind (1st year parts and labour, 2 further years parts only)



OVERALL DIMENSIONS:

57cm – 99cm (H) x 69cm (W) (excluding arms) x 183cm (L)
(maximum in plinth format)

ACCESSORIES:

Paper Roll Holder

Product Code: CHA02

DOHERTY PHLEBOTOMY CHAIR:


PRODUCT CODE: CHE01/(Colour)

DOHERTY PHLEBOTOMY CHAIR

WITH BREATHING HOLE:

PRODUCT CODE: CHE01/BH/(Colour)

DOHERTY VARI-HEIGHT TREATMENT CHAIR



KEY FEATURES:

- Electrically operated variable height range
- Electric handset control for height adjustment
- Breathing hole available as an optional extra
- Split leg sections which are extendable
- The Variable Height Treatment Chair has the added advantage of also being suitable for use as a plinth. Where cost or space constraints limit the amount of equipment that can be purchased this facility makes the chair particularly attractive
- A paper roll holder is available as an "easi fit" accessory
- This treatment chair is fitted with arm pads as standard however Phlebotomy arms can be purchased as an accessory
- Easy lift gas assisted backrest and foot section
- Maximum user weight 160kgs (25 stone)
- Dual sided operation
- Distinctive Doherty style upholstery with rounded corners for hygiene and cleanliness
- Twin retractable castors for ease of movement
- Simple levelling mechanism for uneven floors
- Attractive white epoxy coating with Biocote
- Delivered fully assembled
- Available in 12 modern colour upholstery options
- 3 years warranty for peace of mind (1st year parts and labour ,2 further years parts only)



OVERALL DIMENSIONS:

57cm – 99cm (H) x 69cm (W) (excluding arms) x 184cm (L)
(maximum in plinth format)

ACCESSORIES:

PAPER ROLL HOLDER: CHA02


PHLEBOTOMY ARMS (pair): CHA01

DOHERTY VARIABLE HEIGHT TREATMENT CHAIR:

PRODUCT CODE: CHE02/(Colour)

DOHERTY VARIABLE HEIGHT TREATMENT CHAIR:

PRODUCT CODE: CHE02/BH/(Colour)




DOHERTY FIXED HEIGHT TREATMENT CHAIR



KEY FEATURES:

- This treatment chair is fitted with arm pads as standard, however Phlebotomy arms can be purchased as an accessory
- Split leg sections
- Easy lift gas assisted backrest and foot section
- Maximum user weight 160kgs (25 stone)
- Dual sided operation
- Distinctive Doherty style upholstery with rounded corners for hygiene and cleanliness
- Twin retractable castors for ease of movement
- Simple levelling mechanism for uneven floors
- Attractive white epoxy coating with Biocote
- Delivered fully assembled
- Available in 12 modern colour upholstery options
- 3 years warranty for peace of mind (1st year parts and labour, 2 further years parts only)



OVERALL DIMENSIONS:

60cm (H) x 69cm (W) (excluding arms) x 181.5cm (L)
(maximum with backrest angle at its limit 10 degrees from horizontal)

DOHERTY FIXED HEIGHT TREATMENT CHAIR:

PRODUCT CODE: CHF01/ (Colour)

ACCESSORIES:

PHLEBOTOMY ARMS (pair): CHA01

