

LINEN TROLLEYS



Within the Sidhil range of trolleys is an extensive collection of linen, soiled, changing and cylinder trolleys, which are suitable for hospital, nursing and residential care home use. All trolleys are manufactured to the highest standard and incorporate Biocote® within the epoxy coating to significantly reduce the risk of cross contamination and infection. A wide range of optional accessories and linen bags are available to compliment the range.

LINWOOD SINGLE LINEN TROLLEY (frame only)

- Basic trolley for linen collection and distribution
- 50mm castors
- Optional linen bags

DIMENSIONS:

91.5cm (H) x 48.5cm (Diameter)
46cm x 76cm (Bag size)

Product Code: 1523/2



NYLON BAGS (to fit all Sidhil linen & changing trolleys)

- Available in White, Red, Blue or Green. Add /W, /R, /B or /G to end of product code for colour option
- Product Code: 1523/004/(Letter)



LINCOLN DOUBLE LINEN TROLLEY (frame only)

As above with additional linen bag ring

DIMENSIONS:

91.5cm (H) x 48.5cm (W) 107.5cm (L)
48.5cm (Diameter) Bag ring
46cm x 76cm (Bag size)

Product Code: 1521/2




CHEPSTOW CHANGING TROLLEY (with nylon bag)


- Dual purpose trolley
- Three removable shelves
- 100mm castors
- Supplied with one white linen bag

DIMENSIONS:

104cm (H) x 51cm (W) x 113cm (L)
46cm x 46cm (Tray sizes)
48.5cm (Diameter) Bag ring
46cm x 76cm (Bag size)
23cm (Height to bottom shelf)
91.5cm (Height to top shelf)

Product Code: 1599





CYLINDER TROLLEYS



AFTON CYLINDER TROLLEY

- Suitable for (AF) size cylinders (66cm x 18cm)
- Front release bracket
- Integrated hook
- Anti-static castors

DIMENSIONS:

108cm (H) x 33cm (W) x 34cm (D)
66cm x 18cm (Cylinder size)

Product Code: 1530/I



FALMOUTH CYLINDER TROLLEY

As above, for (F) size cylinders (94cm x 12.5cm)

DIMENSIONS:

108cm (H) x 33cm (W) x 34cm (D)

Product Code: 1527



UNIVERSAL CYLINDER TROLLEY

- Suitable for all size cylinders up to G
- Front release strap
- Integrated hook
- Anti-static castors

DIMENSIONS:

106cm (H) x 30cm (W) x 29cm (D)
Minimum cylinder diameter 9cm
Maximum cylinder diameter 20cm

Product Code: 1531

